

Kodak

SoftWear Lens

for Comfortable Computer Viewing

Avoid Eyestrain with KODAK SoftWear® Lenses

Spending long hours hunched in front of a computer does nothing for posture or eyesight. Computer-based projects or large amounts of detail work demand the use of a task-oriented lens – KODAK SoftWear Lenses.

KODAK SoftWear Lenses are digitally-created, full backside lenses, assisting the eyes for prolonged intermediate and near-viewing tasks. KODAK SoftWear Lenses help maintain productivity in computer-viewing tasks.



Vision First Design™ Technology

KODAK SoftWear Lenses incorporate Vision First Design Technology into a full backside, task-oriented lens. This technology provides KODAK SoftWear Lenses with a smooth power progression and excellent binocular balance to ease patient adaptation and provide excellent visual performance.

Easy To Order and Fit

KODAK SoftWear Lenses are as easy to order and to fit as KODAK Progressive Lenses. With a set computer monitor viewing distance of 24", KODAK SoftWear Lenses simplify lens dispensing without compromising visual performance.

KODAK SoftWear Lenses

Lens Design: Task Oriented/Computer

Viewing Areas: Near and Intermediate

Materials: 1.50, Poly, 1.67, KODAK Total Blue
1.56 and 1.67

Minimum Fitting Height: 17mm



KODAK SoftWear Lenses are digital lenses designed with larger near and intermediate viewing areas.



Greater Levels of Protection

- Add a KODAK Anti-Reflective Lens Coating to reduce lens reflections and eyestrain.
- KODAK Total Blue Lens filters up to 80% of high energy visible blue light allowing your patients to view digital devices safely.

KODAK SoftWear Lenses Availability

Standard Resin

Maximum Cutout: 76mm Rx Range: -8.00 to +4.00

Polycarbonate

Maximum Cutout: 76mm Rx Range: -10.00 to +5.00

1.67 High Index

Maximum Cutout: 76/80mm Rx Range: -12.00 to +10.00

KODAK Total Blue™ 1.56

Maximum Cutout: 80mm Rx Range: -12.00 to +6.00

KODAK Total Blue 1.67

Maximum Cutout: 76/80mm Rx Range: -12.00 to +10.00

Add Powers: +0.75 to +3.50 in .25 D Steps

For Your Patients

Since KODAK SoftWear Lenses devote a majority of the vision support to the near and intermediate viewing areas, please advise your patients that KODAK SoftWear Lenses should not be worn while driving or performing other far-viewing tasks.

See the *Colors* of Life



Signet Armorlite
5803 Newton Drive, Suite A
Carlsbad, CA 92008
800-950-5367
www.signetarmorlite.com



Illustrations for demonstration purposes only.

The Kodak trademark and trade dress are used under license from Kodak by Signet Armorlite, Inc. SoftWear is a registered trademark; Total Blue and Vision First Design are trademarks of Signet Armorlite, Inc.
©2017 Signet Armorlite, Inc.

2/2017 PN 359-196